



FullCare

EXTERIOR VISUALIZATION OF REPIRATOR

The way of connecting between "rubber band" & "mask body" is ultrasonic welding

The size of FFP2 logo is 4.4cm(L) × 1cm(W)

The diameter of connecting point is 0.8cm

Part 1 material: the 1st layer:PP nonwoven fabric 35gsm.
the 2nd layer:FFP2 filter.
the 3rd layer:Pin sub-cotton 240gsm.

The white plastic bar size is:
8.3cm(L)×0.8cm(W)×0.05cm(T)

PART 1

PART 2

The front yellow rubber size is 30cm(L)×
0.8cm(W)×0.04cm(T)

PART 4

PART 3

REVISION 1

DRAWING No. :

PART 1	Respirator Body	
PART 2	Rubber Bands	The elastic modulus of the straps is 0.7 N/mm.
PART 3	Nose Bar	
PART 4	Logo	

Mark	Content	Person	Date	Audit	Weight	Scale
	Design				35g	1:1
1	Proofread				2/4 Page	

